Leukomed® Sorbact®

EN

Instructions for Use

Read the instructions for use before starting treatment with Leukomed® Sorbact®.

Device Description

Leukomed® Sorbact® is a hydrophobic microbe-binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACC™-coated wound contact layer combined with a wound pad and a transparent film coated with acrylic adhesive. Leukomed® Sorbact® absorbs and retains exudate and enables a moist wound environment. The waterproof backing film provides protection against external contamination and allows excess fluid to evaporate. Leukomed® Sorbact® manages wound infections.

Mode of Action

Sorbact® can bind hydrophobic microbes based on in vitro testing of the following species: *S. aureus* including MRSA, *S. gordonii*, *P. aeruginosa*, and *E. coli*. When the dressing is removed, microbes adhered to the dressing will be removed.

Intended Use

Leukomed® Sorbact® is indicated for use in the management of dry to low exuding wounds, as well as partial to full thickness wounds (including clean, contaminated, colonized or infected wounds): post-operative wounds and trauma wounds.

Contraindications

There are no known contraindications for the use of Leukomed® Sorbact®. Warnings and precautions should be noted.

Warnings and Precautions

Do not re-use. Leukomed® Sorbact® is for single patient and single use only. Reusage may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate, polyurethane, polyolefin, polyacrylate, polyester and viscose). In common with all adhesive products, cases of irritation to the surrounding skin may occur. It should also be noted that inappropriate use or too frequent dressing changes, particularly in patients with fragile skin, may result in superficial

skin trauma. Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of hydrophobic microbes. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use

1. Prepare the wound and surrounding skin according to local clinical practice. Ensure that the surrounding skin is clean and dry.
2. Select an appropriate dressing size for the wound. Do not cut the dressing.
3. Remove the dressing from the peel pouch using an aseptic technique.
4. Remove the adhesive protective liners marked (1) from the wound contact side of the dressing and apply the dressing over the wound. Ensure that the Sorbact® wound contact layer comes into direct contact with the complete wound surface.
5. Press the borders of the dressing to the surrounding skin and remove the stabilizing film marked (2) on top. Start at the green tab.
6. The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
7. For easy removal, it is recommended to remove the dressing by pulling and stretching the film of the dressing parallel to the skin. Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound. Leukomed® Sorbact® is MRI safe. Prior to radiation therapy, remove Leukomed® Sorbact® if placed in the radiation field. A new dressing can be applied following treatment.

Storage and Disposal

Leukomed® Sorbact® shall be stored in a dry place and kept away from sunlight. Disposal should be made according to local environmental procedures.

Caution

Federal law (U.S.A.) restricts the sale of this device except by a physician's orders.

Leukomed® and Leukoplast® are registered trademarks of BSN medical GmbH. Sorbact® is a registered trademark of ABIGO Medical AB. Leukomed®, Leukoplast® son marcas comercial registradas de BSN medical GmbH. Sorbact® es una marca comercial registrada de ABIGO Medical AB.

Made in Germany

® = registered trademark

Patent information: /
www.essity.com/patents

Rx only

Not made with natural rubber latex
No fabricado con látex de caucho natural

REF

Catalogue number

LOT

Batch code

Use-by date

Date of manufacture

UDI

Unique device identifier

MD

Medical device

Caution

Consult instructions for use

Do not re-use

Keep away from sunlight

Keep dry



Do not use if package is damaged and consult instructions for use

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