Instructions for Use

en

Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel

Read the instructions for use before starting treatment with Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel.

Device Description

Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is a bacteria and fungi binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACC™-coated wound contact layer combined with an absorbent polyurethane foam that contains superabsorbent stripes. The product absorbs and retains exudate. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel prevents and treats wound infections and facilitates the wound healing process. Sorbact® wound contact layer is shown to bind wound related pathogens up to seven days *in vitro*.

Intended Purpose

Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is intended for use in management of clean, contaminated, colonized or infected low to moderate exuding wounds, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic foot ulcers and leg ulcers. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is intended to be used on superficial wounds.

Intended User and Use Environment

Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is intended to be used on children and adults. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is intended to be used by healthcare professionals, or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

Mode of Action

Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

Clinical Benefits

Use of Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel can reduce bioburden and prevent infection. The dressing helps reducing pain and can improve wound healing and decrease wound size. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel absorbs and retains exudate. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel reduces the risk of maceration and enables a moist wound environment.

Contraindications

There are no known contraindications to the use of Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel. Warnings and precautions should be noted.

Warnings and Precautions

Do not re-use. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate, polyurethane, polyolefin and sodium polyacrylate). Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of microorganisms. Do not use with oxidizing solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent component of the dressing. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use

- Prepare the wound and surrounding skin according to local clinical practice. Ensure that the surrounding skin is clean and dry.
- Choose an appropriate dressing size that covers the wound and overlaps the wound edges. **(A)**
- Remove the dressing from the peel pouch using an aseptic technique.
- If required, the dressing may be cut using an aseptic technique to suit various wound shapes and locations. Discard unused pieces of the dressing.
- Ensure that the Sorbact® wound contact layer comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing. **(B)**
- Fixate with an elastic bandage or skin friendly tape.
- The dressing changing frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
- Gently remove any additional fixation and the dressing from the wound and discard. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel can be used in conjunction with compression therapy. Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel is MRI safe. Prior to radiation therapy, remove Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel if it is placed in the radiation field. A new dressing can be applied following treatment.

Storage and Disposal

Cutimed® Siltec® Sorbact® / Cutimed® Siltec® Sorbact® Heel shall be stored dry and kept away from sunlight. Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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