EN

Read the instructions for use before starting treatment with Sorbact® Gel Compress.

Device Description

Sorbact® Gel Compress is a non-absorbent bacteria and fungi binding wound dressing, based on Sorbact Technology. It consists of a Sorbact® DACC™ -coated wound contact layer covered with a hydrogel containing carbomer and propylene glycol (10%). Sorbact® Gel Compress treats and helps to prevent wound infections and facilitates the wound healing process. Sorbact® wound contact layer is shown to bind wound related pathogens up to seven days in vitro.

Intended Purpose

Sorbact® Gel Compress is intended for use in management of clean, contaminated, colonized or infected dry to exuding wounds, such as surgical wounds, traumatic wounds, burns, pressure ulcers, diabetic foot ulcers and leg ulcers. Sorbact® Gel Compress can be used on both superficial and deep wounds.

Intended User and Use Environment

Sorbact® Gel Compress is intended to be used on children and adults. Sorbact® Gel Compress is intended to be used by healthcare professionals, or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

Mode of Action

Sorbact® Gel Compress binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown in vitro. There is no release of antimicrobial agents to the wound.

Clinical Benefits

Use of Sorbact® Gel Compress can reduce bioburden and helps to prevent infection. The dressing helps reducing pain and can improve wound healing and decrease wound size. Sorbact® Gel Compress donates moisture, enables a moist environment and facilitates autolytic debridement.

Contraindications

There are no known contraindications to the use of Sorbact® Gel Compress. Warning and precautions should be noted.

Warnings and Precautions

Do not re-use. Sorbact® Gel Compress is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection.

Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed.

Do not use on patients with a known allergy to propylene glycol or with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate, carbomer, propylene glycol and water).

Caution should be taken when exposing children below 5 years and patients suffering from liver and kidney disease, due to limited capacity to metabolize propylene glycol.

High exposure, e.g. in case of simultaneous use of other products containing propylene glycol, may in rare cases lead to adverse events, e.g. hypersmolarity or metabolic acidosis.

Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Do not re-sterilize.

If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use

1. Prepare the wound and surrounding skin according to local clinical practice. **2.** Select an appropriate dressing size. **3.** Remove the dressing from the peel pouch using an aseptic technique. **4.** If the dressing is cut, use an aseptic technique. **5.** Apply the dressing. Ensure that the dressing comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing. The dressing should not overlap the surrounding skin. Avoid dense packing when used in deep wounds. **6.** Apply a moisture preserving secondary dressing appropriate for the exudate level, such as foam dressing, absorbent pad or film dressing, and fixate. **7.** The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days. **8.** Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound. Sorbact® Gel Compress can be used in conjunction with compression therapy. Sorbact® Gel Compress is MRI safe. Prior to radiation therapy, remove Sorbact® Gel Compress if placed in the radiation field. A new dressing can be applied following treatment.

Storage and Disposal

Sorbact® Gel Compress shall be stored dry and kept away from sunlight. Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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