Cutimed® Sorbact Ribbon®

(en)

Read the instructions for use before starting treatment with Cutimed® Sorbact Ribbon®.

Device Description

Cutimed® Sorbact Ribbon® is a DACC™-coated bacteria and fungi binding ribbon gauze, based on Sorbact® Technology. Cutimed® Sorbact Ribbon® prevents and treats wound infections and facilitates the wound healing process. It is suitable as a wound filler and allows passage of wound exudate into a secondary dressing. Cutimed® Sorbact Ribbon® can also be placed in skin folds to treat fungal infections (intertrigo).

Intended Purpose

Cutimed® Sorbact Ribbon® is intended for use in management of clean, contaminated, colonized or infected exuding wounds, such as surgical wounds, traumatic wounds, pressure ulcers, diabetic foot ulcers and leg ulcers. Cutimed® Sorbact Ribbon® can be used on superficial and deep wounds such as cavity wounds and fistulas. Cutimed® Sorbact Ribbon® is also used to treat fungal infections in skin folds (intertrigo).

Intended User and Use Environment

Cutimed® Sorbact Ribbon® is intended to be used on children and adults. Cutimed® Sorbact Ribbon® is intended to be used by healthcare professionals or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

Mode of Action

Cutimed® Sorbact Ribbon® binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

Clinical Benefits

Use of Cutimed® Sorbact Ribbon® can reduce bioburden and prevent infection. The dressing helps reducing odor and pain, and can improve wound healing and decrease wound size. Cutimed® Sorbact Ribbon® can also treat fungal infections in skin folds (intertrigo).

Contraindications

There are no known contraindications to the use of Cutimed® Sorbact Ribbon®. Warnings and precautions should be noted.

Warnings and Precautions

Do not re-use, Cutimed® Sorbact Ribbon® is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cotton). Do not use in combination with fatty products, such as ointments, creams and solutions, including topical antifungal treatments, as they may decrease the binding of microorganisms. Do not wrap around extremities, as circulation may be compromised. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

Instructions for Use for management of wounds

- Prepare the wound and surrounding skin according to local clinical practice.
- Select an appropriate dressing size for the wound. The dressing may overlap the wound margins if needed.
- Remove the dressing from the peel pouch using an aseptic technique.
- If needed, cut across the dressing using an aseptic technique. Do not cut over the wound. Discard any unused dressing.
- Apply the dressing, ensuring that the dressing comes into direct contact with the complete wound surface to allow microorganisms to bind to the dressing. In cavities or fistulas, avoid dense packing. Always leave a part of the dressing outside the cavity/fistula for easy removal. If the dressing is cut, always leave the cut end outside the wound to avoid residual threads.
- Apply a secondary dressing appropriate for the exudate level such as a foam dressing, absorbent pad or gelling fiber dressing, and fixate.
- The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
- Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of wound bed and surrounding tissue.

Cutimed® Sorbact Ribbon® can be used in conjunction with compression therapy. Cutimed® Sorbact Ribbon® is MRI safe. Cutimed® Sorbact Ribbon® does not contain x-ray threat. Prior to radiation therapy, remove Cutimed® Sorbact Ribbon® if placed in the radiation field. A new dressing can be applied following treatment.

Instructions for Use for treatment of fungal infections in skin folds

- Prepare the affected area according to local clinical practice.
- Select an appropriate dressing size for the affected area. If needed cut across the dressing to a suitable size.
- Apply the dressing. Ensure that the dressing comes into direct contact with the affected area to allow microorganisms to bind to the dressing.
- For fixation, undervear or a skin friendly adhesive can be used. The adhesive should be applied outside the affected area.
- Daily dressing change in conjunction with inspection and regular hygiene routines is recommended.
- Continue to use the dressing for 2–3 days after the symptoms have resolved to avoid recurrence.

Storage and Disposal

Cutimed® Sorbact Ribbon® shall be stored dry and kept away from sunlight.

Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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